myoPRESSURE DPM/DPM-T Quick Start

STEP 1: Home/Start Screen:

- 1 Select the module myoPRESSURE
- 2 Select/create a subject
 - Enter Subject Name and other subject properties
- 3 Select Edit or New Configuration
 - Edit will re-enter an already existing configuration.
 - New will allow you to create a new configuration.
 - If you mark a configuration from the list and hit the green
 Measure button you will directly enter the Measure Menu.



STEP 2: Create or edit a measurement configuration

- Insert DPM-T (treadmill) or DPM (floor plate)
 - Drag the device icon from the right tool bar to the configuration window.

Configure the platform

- For DPM/DPM-T select Treadmill mode for regular gait analysis and Stance mode for unilateral or bilateral stance analysis.
- Unless you own high speed plates, the recommended sampling rate is 100 Hz.

Continue with

2


▶ Start Measure



STEP 3: Calibrate and Measure

- 1 Calibrate the platform
 - Click on Continue to calibrate the unloaded platform (subject is not standing on platform).
 - Let the subject stand/walk on the platform to check if real time pressure prints are
- 2 Click on Record to start the measurement

Continue Please step off the platform and press 'Continue'.



- 3 Perform the test activity and press 'Stop' to finish the recording.
 - | Sea-Date | Sea |

- 3 Click on "Stop" to stop the recording
- 4 Enter a record name and save the record
 - Save & View will save the record and move it to the record viewer for viewing
 - Save & measure will immediately go back to the measurement window.
 Discard & measure will delete the record and return to the measurement

myoPRESSURE DPM/DPM-T Quick Start

STEP 4: Viewer

Replay the record 1

Use video play buttons below the time line

Review channels

- Double click on channel names to display them in separate windows, or drag the channel by its name tab to the plus sign.
- Mouse drag a channel over another channel to overlay the two (or more) signals.

Use layout system to arrange signals on screen

· Store and load screen configurations

3

Continue with





STEP 5: Report selection, channel selection and analysis period definition

Select a report in tab myoPRESSURE Reports

Preview/study report contents and report definition

- Each report comes with a pre-configured mode for analysis period definition, which is
- Study the details of a chosen period definition, which may be done with marker sets. foot step detection, or min/max angle event detection algorithms.

Click on OK to enter the Report channel selection sub step

Select channels

- Analysis channels can be selected via Group, Type, Body side, Sensor name, or physical channel number. For DPM-T/DPM analysis the 2 force channels must be selected.
- This channel filter is especially helpful if pressure recording is combined with other modules/devices like EMG, 3D motion, and/or video.

Select analysis period 4

- myoPressure DPM-T/DPM works with an automatic step and analysis period detection.
- Use the mouse to place markers that limit analysis to a special portion within the record and click Set in right tool bar.
- The colored bars in time line indicate the gait phases (red & green colored).

Click on Next to continue to go to the analysis report





STEP 6: Read and print/export a report

Scroll through the report results

Use Print, View, or Analysis Options

- Reports can be printed, copied to clipboard, or exported to several files formats.
- The report can be reanalyzed to adjust channel selection or the period definition(s).
- The report can be compared to another record/report of the same record format. Averaged signals can be stored as time normalized curves in normative or group data.
- Click on Home to go back to home screen or click Measure in the

Top navigation bar to directly measure with the same measurement configuration

