

Timed Up & Go Report



Patient
Project
Last Name
First Name


2023 Demo Records
TUG
Demo Records

Record
Name
Date Measured
Number of periods

TUG (Limp)
1/26/2023 16:21
1



Protocol



Results

Timing, s	Total	Stand	Walk 1	Turn 1	Walk 2	Turn 2	Sit
	15.6	1.5	4.7	2.0	5.1	0.9	1.4
Trunk Flex. ROM, deg				41.5			
Trunk Lat. ROM, deg				3.6			
Walk Cadence, steps/min				101			
Walk 1st Double Support, %	Left			16.2			
	Right			12.4			
Walk Stance Ratio (R - L) / (R + L), %				-1.6			
Walk Foot Pitch At Foot Off, deg	Left			19.9			
	Right			38.3			
Turn 1 Step Count, steps				2			
Turn 1 Trunk Rotation Speed Max, deg/s				24.7			

Subject Comments

Report Comments